

SPECIFICATION

Strontium Nitrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay[Sr(NO ₃) ₂] Solubility in water pH of a 50g/L solution at 25°C Chloride (Cl) Sulfate (SO ₄) Sodium (Na) Potassium (K) Magnesium (Mg) Calcium (Ca) Barium (Ba) Lead (Pb) Iron (Fe)	White crystalline powder 98.0~102.0%(mass/mass) to pass test 5.0~7.0 max.0.001%(mass/mass) max.0.003%(mass/mass) max.0.01%(mass/mass) max.0.01%(mass/mass) max.0.005%(mass/mass) max.0.02%(mass/mass) max.0.5%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass)

JIS K 8554:2020(*:Additional test performed by Wako but not required by JIS)