

## SPECIFICATION

Disodium Hydrogenphosphate  
for pH Standard Solution

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (Na <sub>2</sub> HPO <sub>4</sub> ) (after drying at 110°C)	min.99.5%(mass/mass)
Solubility in water	to pass test
Loss on drying at 110°C	max.0.5%(mass/mass)
Chloride (Cl)	max.0.002%(mass/mass)
Nitrate	to pass test
Condensed phosphate	to pass test
Sulfate (SO <sub>4</sub> )	max.0.005%(mass/mass)
Potassium (K)	max.0.04%(mass/mass)
Calcium (Ca)	max.0.03%(mass/mass)
Lead (Pb)	max.0.002%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.002%(mass/mass)
Ammonium (NH <sub>4</sub> )	max.0.001%(mass/mass)

JIS K 9020:2021(\*:Additional test performed by Wako but not required by JIS)