

Code No.19-0728

## SPECIFICATION

pH5.10 Sodium Citrate Buffer Solution (5 times concentrated)  
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* pH at 25°C	5.05~5.15
Concentration of Cl (Cl <sup>-</sup> )	0.40~0.80mol/L
Concentration of Na (Na <sup>+</sup> )	2.10~2.40mol/L
* Ninhydrin reactive substances	to pass test

\*:Specification of 5 times diluted solution