

Code No.19-0749

## SPECIFICATION

### Disodium Hydrogenphosphate 12-Water for Oxidant Analysis for Oxidant Analysis

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay	99.0~100.5%
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	9.0~9.4
Nitrogen compounds (as N)	max.0.001%
Heavy metals (as Pb)	max.5ppm
Suitability for oxidant determination	to pass test

JIS K 8828:2004

(\*:Additional test performed by Wako but not required by JIS)