

Code No.19-1007

SPECIFICATION

Siduron Standard (mixture of isomers)

for Pesticide Residue Analysis

| REQUIREMENT | SPECIFICATION |
|--|----------------------------------|
| Appearance | White, crystals ~ powder or mass |
| Solubility in methanol | to pass test |
| Melting point | 136~144℃ |
| Loss on drying over silica gel in vacuum | max.0.1% |
| Assay (HPLC) (mixture of isomers) | min.98.0% |
| Assay (qNMR) (mixture of isomers) | min.98.0% |