

Code No.19-1099

SPECIFICATION

Stearic Acid
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, flakes or crystalline powder ~ powder
Assay (C18H36O2) (GC)	min.95.0%(mass/mass)
Melting point	68~71℃
Residue after ignition (as sulfate)	max.0.01%(mass/mass)
Acid value	195~200
Iodine value	max.1.5

JIS K 8585:2015(*:Additional test performed by Wako but not required by JIS)