

## SPECIFICATION

## Sodium Azide

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ nearly white, crystalline powder or mass min.98.0%(mass/mass)  to pass test  max.0.2%(mass/mass)  max.0.1%(mass/mass)  max.2ppm(mass/mass)
Assay (NaN <sub>3</sub> )	
Solubility in water	
Alkalinity (as NaOH)	
Loss on drying at 105°C	
Chromium (Cr)	

JIS K 9501 : 2019(\*:Additional test performed by Wako but not required by JIS)