

## SPECIFICATION

Sample Buffer Solution (2ME-)(X4)  
for Electrophoresis

| REQUIREMENT                       | SPECIFICATION     |
|-----------------------------------|-------------------|
| Appearance                        | Blue clear liquid |
| Density (20°C)                    | 1.109~1.115g/mL   |
| Absorption maximum                | 590~593nm         |
| Absorbance (3→100, $\lambda$ max) | 0.55~0.75         |
| pH at 25°C                        | 6.7~6.9           |
| Practical test (SDS-PAGE)         | to pass test      |