

## SPECIFICATION

## Sodium L-Tartrate

## Japanese Pharmaceutical Excipients

| REQUIREMENT  | SPECIFICATION  |
|--|--|
| * Description<br>Identification<br>Optical rotation [ $\alpha$ ]D20<br>pH<br>Purity<br>(1) Clarity and color of solution<br>(2) Sulfate<br>(3) Oxalate<br>(4) Heavy metals<br>(5) Arsenic<br>Loss on drying<br>Assay (as anhydrous) (after drying) | Colorless crystals or white crystalline powder<br><br>to pass test<br>+25.0~+27.5°<br>7.0~9.0<br>-<br>to pass test<br>not more than 0.019%<br>to pass test<br>not more than 10ppm<br>not more than 2ppm<br>14.0~17.0%<br>not less than 98.5% |

Japanese Pharmaceutical Excipients(\*:Reference item)