

Code No.19-1508

## SPECIFICATION

### Sarcosine Standard

for Amino Acid Analysis

| REQUIREMENT                               | SPECIFICATION                                                       |
|-------------------------------------------|---------------------------------------------------------------------|
| Appearance                                | White ~ nearly white,<br>crystalline powder or mass<br>to pass test |
| Solubility in water                       |                                                                     |
| Loss on drying at 105°C                   | under investigation                                                 |
| Suitability for amino acids determination | under investigation                                                 |
| Assay (volumetric analysis)               | min.98.0%                                                           |
| Assay (subtracting method)                | under investigation                                                 |