

Code No.19-1617

SPECIFICATION

Saikosaponin b2

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and TLC)

| REQUIREMENT | SPECIFICATION |
|--|--------------------------------------|
| Description (JP method) | White crystals or crystalline powder |
| Identification (JP method) | to pass test |
| Purity (JP method) | - |
| Related substances (JP method) | to pass test |
| Unity of peak (JP method) | to pass test |
| Assay (qNMR) (JP method) | under investigation |
| * Solubility in methanol | to pass test |
| * Loss on drying over silica gel in vacuum | max.10.0% |

This product is suited for Saikosaponin b2 for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (*:Additional test performed by Wako)