

Code No.20-0008

## SPECIFICATION

## DL-Tartaric Acid

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystalline powder ~ powder
Solubility in water	to pass test
Solubility in ethanol	to pass test
Specific rotation [ $\alpha$ ]D20(c=20,H2O)	-0.2~+0.2°
Residue after ignition (as sulfate)	max.0.1%
Chloride (Cl)	max.0.002%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	min.99.0%