

SPECIFICATION

Tin(II) Chloride Dihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (SnCl ₂ ·2H ₂ O)	min.97.0%(mass/mass)
Solubility in hydrochloric acid	to pass test
Sulfate (SO ₄)	max.0.002%(mass/mass)
Sodium (Na)	max.0.005%(mass/mass)
Potassium (K)	max.0.002%(mass/mass)
Copper (Cu)	max.0.001%(mass/mass)
Calcium (Ca)	max.0.002%(mass/mass)
Lead (Pb)	max.0.003%(mass/mass)
Arsenic (As)	max.2ppm(mass/mass)
Iron (Fe)	max.0.002%(mass/mass)
Ammonium (NH ₄)	max.0.002%(mass/mass)

JIS K 8136:2017*:Additional test performed by Wako but not required by JIS