

Code No.20-0187

## SPECIFICATION

## Toluene

Wako 1st Grade

| REQUIREMENT                          | SPECIFICATION          |
|--------------------------------------|------------------------|
| Appearance                           | Colorless clear liquid |
| Density (20°C)                       | 0.864~0.868g/mL        |
| Refractive index 20°C                | 1.495~1.498            |
| Water                                | max.0.05%              |
| Sulfur compounds                     | to pass test           |
| Thiophene                            | to pass test           |
| Substances darkened by sulfuric acid | to pass test           |
| Assay (capillary GC)                 | min.99.0%              |