

## SPECIFICATION

## Trichloroacetic Acid

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CCl <sub>3</sub> COOH) (after drying) Solubility in water Loss on drying Residue after ignition (as sulfate) Chloride (Cl) Nitrate Phosphate (PO <sub>4</sub> ) Sulfate (SO <sub>4</sub> ) Lead (Pb) Iron (Fe) Substances darkened by sulfuric acid	White, crystals ~ crystalline powder or mass min.99.0%(mass/mass) to pass test max.1.0%(mass/mass) max.0.02%(mass/mass) max.0.001%(mass/mass) to pass test max.5ppm(mass/mass) max.0.01%(mass/mass) max.0.001%(mass/mass) max.5ppm(mass/mass) to pass test

JIS K 8667 : 2019(\*:Additional test performed by Wako but not required by JIS)