

Code No.20-0288

## SPECIFICATION

## 2,2,4-Trimethylpentane

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay [CH <sub>3</sub> C(CH <sub>3</sub> ) <sub>2</sub> CH <sub>2</sub> CH(CH <sub>3</sub> )CH <sub>3</sub> ] (GC)	min.99.0%(mass/mass)
Density (20°C)	0.690~0.693g/ml
Refractive index 20°C	1.390~1.393
Residue after evaporation	max.0.001%(mass/mass)
Substances darkened by sulfuric acid	to pass test

JIS K 9703:2013(\*:Additional test performed by Wako but not required by JIS)