

Code No.20-0706

## SPECIFICATION

## 1,1,2,2-Tetrachloroethane

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless ~ nearly colorless, clear liquid
Assay (Cl <sub>2</sub> CHCHCl <sub>2</sub> )(GC)	min.97.0%(area/area)
Density (20°C)	1.593~1.601g/mL
Refractive index 20°C	1.493~1.497
Water	max.0.1%(mass/mass)
Acidity (as HCl)	max.0.003%(mass/mass)

JIS K 9520:2024(\*:Additional test performed by Wako but not required by JIS)