

Code No.20-0874

SPECIFICATION

Tetrahydrofuran
Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|------------------------|
| * Appearance | Colorless clear liquid |
| Assay (C ₄ H ₈ O) (GC) | min.99.5%(area/area) |
| Density (20℃) | 0.884~0.889g/mL |
| Refractive index 20℃ | 1.406~1.409 |
| Water | max.0.05%(mass/mass) |
| Peroxides (as H ₂ O ₂) | max.0.005%(mass/mass) |

JIS K 9705:2024*:Additional test performed by Wako but not required by JIS(with antioxidants)