

## SPECIFICATION

0.5mol/L Tetrabutylammonium Dihydrogenphosphate Solution  
for Ion-pair Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (200nm)	max.0.50
Absorbance (210nm)	max.0.30
Absorbance (240nm)	max.0.15
Absorbance (254nm)	max.0.15
Absorbance (300nm)	max.0.10
Absorbance (350nm)	max.0.05
pH at 25°C	under investigation
Fluorescence test	to pass test
Concentration	0.45~0.55mol/L