

Code No.20-1694

## SPECIFICATION

0.25w/v% Trypsin-1mmol/L EDTA · 4Na Solution with Phenol Red  
for Cell Culture

| REQUIREMENT     | SPECIFICATION    |
|-----------------|------------------|
| Appearance      | Red clear liquid |
| pH              | 7.0~7.6          |
| Osmolality      | 275~325mOsm/kg   |
| Sterility test  | to pass test     |
| Mycoplasma test | to pass test     |
| Practical test  | to pass test     |