

SPECIFICATION

Thymol for spray solution for the Japanese Pharmacopoeia General Tests

| REQUIREMENT | SPECIFICATION |
|----------------------------|--------------------------------------|
| Description (JP method) | White crystals or crystalline powder |
| Identification (JP method) | to pass test |
| Melting point (JP method) | 49~52°C |
| Purity (JP method) | - |
| Foreign phenol compounds | to pass test |
| Solubility in ethanol | to pass test |
| Practical test | to pass test |
| Assay (capillary GC) | min.98.0% |