

## SPECIFICATION

 $(\pm)$ - $\alpha$ -Tocopherol Acetate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ yellow, clear liquid
Refractive index 20°C	1.494~1.499
Specific absorbance E1%1cm (284nm)	41.0~45.0
Water	max.0.2%
Acidity (as CH <sub>3</sub> COOH)	max.0.2%
$\alpha$ -Tocopherol	max.0.5%
Assay (HPLC)	min.96.0%