

Code No.20-1870

## SPECIFICATION

### Tetrahydrofuran, Deoxidized, with Stabilizer for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.884~0.889g/mL
Refractive index 20°C	1.406~1.409
Water	max.0.001%
Dissolved oxygen	max.1ppm
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.005%
Assay (capillary GC)	min.99.5%