

Code No.20-1933

## SPECIFICATION

## Tylosin Standard

for High Performance Liquid Chromatography

| REQUIREMENT                                  | SPECIFICATION                                 |
|--|---|
| Appearance                                   | White ~ slightly yellowish brown, crystalline |
| Solubility in acetonitrile                   | powder ~ powder<br>to pass test               |
| Specific rotation [ $\alpha$ ]D20(c=1,CH3OH) | under investigation                           |
| Loss on drying at 105°C                      | under investigation                           |
| Assay (HPLC)                                 | min.97.0%                                     |