

## SPECIFICATION

Tetrahydrofuran, with Stabilizer  
for HPLC

| REQUIREMENT                                  | SPECIFICATION          |
|----------------------------------------------|------------------------|
| Appearance                                   | Colorless clear liquid |
| Density (20°C)                               | 0.884~0.889g/mL        |
| Refractive index 20°C                        | 1.406~1.409            |
| Absorbance (300nm)                           | max.1.00               |
| Absorbance (320~600nm)                       | max.0.10               |
| Water                                        | max.0.05%              |
| Peroxide (as H <sub>2</sub> O <sub>2</sub> ) | max.0.005%             |
| Assay (capillary GC)                         | min.99.5%              |
| Stabilizer (BHT)                             | 0.01~0.03%             |