

Code No.21-0009

SPECIFICATION

Fluorescein Sodium

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--|
| * Appearance | Green or yellowish red ~ reddish brown, crystalline powder ~ powder or mass to pass test |
| Solubility in water | |
| Ethanol-insoluble matter | max.1.0%(mass/mass) |
| Absorbance (3mg/L, pH7.0) | min.0.45 |
| Loss on drying at 105°C | max.7.0%(mass/mass) |
| Residue after ignition (as sulfate) (after drying) | 35.0~40.5%(mass/mass) |
| Chloride (Cl) | max.0.1%(mass/mass) |
| Zinc (Zn) | max.0.05%(mass/mass) |
| Sensitivity | to pass test |

JIS K 8830 : 2023(*:Additional test performed by Wako but not required by JIS)