

## SPECIFICATION

Wakogel(R) C-300E  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White powder
Particle size (45~75 $\mu$ m)	min.75%
Particle size (coarser than 75 $\mu$ m)	max.5%
Precipitation volume	2.0~3.0ml/g
Loss on drying at 180°C	max.6.0%
pH (50g/l aqueous slurry,25°C)	5.5~7.5
Flowing rate	0.5~1.4ml/min
Separation power	to pass test
Activity	to pass test