

Code No.24-0046

## SPECIFICATION

## Xylene Cyanol FF

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Water-insoluble matter Ethanol-insoluble matter Absorbance (5mg/l) Loss on drying at 105°C Sensitivity	Dark green ~ dark purple, crystalline powder max.0.5%  max.1.0%  min.0.50  max.5.0%  to pass test

JIS K 8272-1994

(\*:Additional test performed by Wako but not required by JIS)