

Code No.24-0046

## SPECIFICATION

## Xylene Cyanol FF

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Dark green ~ dark purple, crystalline powder
Water-insoluble matter	max.0.5%
Ethanol-insoluble matter	max.1.0%
Absorbance (5mg/l)	min.0.50
Loss on drying at 105°C	max.5.0%
Sensitivity	to pass test

JIS K 8272-1994

(\*:Additional test performed by Wako but not required by JIS)