

SPECIFICATION

Zinc

Reference Material for Volumetric Analysis

REQUIREMENT	SPECIFICATION
* Appearance Assay Copper (Cu) Cadmium (Cd) Tin (Sn) Lead (Pb) Iron (Fe)	Bluish white turnings min.99.99%(mass/mass) max.0.0010%(mass/mass) max.0.0010%(mass/mass) max.5ppm(mass/mass) max.0.0020%(mass/mass) max.0.0010%(mass/mass)

JIS K 8005:2006

*:Additional test performed by Wako but not required by JIS