

Cloud-Clone社「新型コロナウイルス感染研究ターゲットELISA KIT」30%FFキャンペーン



Cloud-Clone Corp.

対象製品価格表 (2020年 10/1~12/25)

製品コード	品名	容量	希望納入価格 (円)	キャンペーン価格 (円)	ターゲット 因子	測定動物種
SEA049Bo	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	175,800	123,000	IFN γ	Bovine
SEA049Ca	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	167,400	117,100	IFN γ	Canine
SEA049Eq	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	180,000	126,000	IFN γ	Equine
SEA049Ga	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	159,100	111,300	IFN γ	Chicken
SEA049Gu	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	167,400	117,100	IFN γ	Guinea pig
SEA049Hu	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	113,500	79,400	IFN γ	Human
SEA049Mu	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	97,200	68,000	IFN γ	Mouse
SEA049Po	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	115,800	81,000	IFN γ	Porcine
SEA049Ra	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	108,200	75,700	IFN γ	Rat
SEA049Rb	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	159,100	111,300	IFN γ	Rabbit
SEA049Si	ELISA Kit for Interferon Gamma (IFNγ)	96Tests	188,400	131,800	IFN γ	Rhesus monkey
SEA056Bo	ELISA Kit for Interleukin 10 (IL10)	96Tests	175,800	123,000	IL10	Bovine
SEA056Ca	ELISA Kit for Interleukin 10 (IL10)	96Tests	167,400	117,100	IL10	Canine
SEA056Cp	ELISA Kit for Interleukin 10 (IL10)	96Tests	180,000	126,000	IL10	Caprine
SEA056Eq	ELISA Kit for Interleukin 10 (IL10)	96Tests	180,000	126,000	IL10	Equine
SEA056Ga	ELISA Kit for Interleukin 10 (IL10)	96Tests	159,100	111,300	IL10	Chicken
SEA056Gu	ELISA Kit for Interleukin 10 (IL10)	96Tests	167,400	117,100	IL10	Guinea pig
SEA056Hu	ELISA Kit for Interleukin 10 (IL10)	96Tests	114,000	79,800	IL10	Human
SEA056Mu	ELISA Kit for Interleukin 10 (IL10)	96Tests	117,200	82,000	IL10	Mouse
SEA056Po	ELISA Kit for Interleukin 10 (IL10)	96Tests	134,400	94,000	IL10	Porcine
SEA056Ra	ELISA Kit for Interleukin 10 (IL10)	96Tests	108,800	76,100	IL10	Rat
SEA056Rb	ELISA Kit for Interleukin 10 (IL10)	96Tests	159,100	111,300	IL10	Rabbit
SEA056Si	ELISA Kit for Interleukin 10 (IL10)	96Tests	188,400	131,800	IL10	Rhesus monkey
SEA061Bo	ELISA Kit for Interleukin 15 (IL15)	96Tests	175,800	123,000	IL15	Bovine
SEA061Hu	ELISA Kit for Interleukin 15 (IL15)	96Tests	108,900	76,200	IL15	Human
SEA061Mu	ELISA Kit for Interleukin 15 (IL15)	96Tests	150,700	105,400	IL15	Mouse
SEA061Ra	ELISA Kit for Interleukin 15 (IL15)	96Tests	159,100	111,300	IL15	Rat
SEA063Bo	ELISA Kit for Interleukin 17 (IL17)	96Tests	175,800	123,000	IL17	Bovine
SEA063Ca	ELISA Kit for Interleukin 17 (IL17)	96Tests	167,400	117,100	IL17	Canine
SEA063Eq	ELISA Kit for Interleukin 17 (IL17)	96Tests	180,000	126,000	IL17	Equine
SEA063Hu	ELISA Kit for Interleukin 17 (IL17)	96Tests	146,500	102,500	IL17	Human
SEA063Mu	ELISA Kit for Interleukin 17 (IL17)	96Tests	108,900	76,200	IL17	Mouse
SEA063Ov	ELISA Kit for Interleukin 17 (IL17)	96Tests	188,400	131,800	IL17	Ovine
SEA063Po	ELISA Kit for Interleukin 17 (IL17)	96Tests	175,800	123,000	IL17	Porcine
SEA063Ra	ELISA Kit for Interleukin 17 (IL17)	96Tests	159,100	111,300	IL17	Rat
SEA079Bo	ELISA Kit for Interleukin 6 (IL6)	96Tests	175,800	123,000	IL6	Bovine
SEA079Ca	ELISA Kit for Interleukin 6 (IL6)	96Tests	167,400	117,100	IL6	Canine
SEA079Cp	ELISA Kit for Interleukin 6 (IL6)	96Tests	160,000	112,000	IL6	Caprine
SEA079Eq	ELISA Kit for Interleukin 6 (IL6)	96Tests	180,000	126,000	IL6	Equine
SEA079Ga	ELISA Kit for Interleukin 6 (IL6)	96Tests	159,100	111,300	IL6	Chicken
SEA079Gu	ELISA Kit for Interleukin 6 (IL6)	96Tests	167,400	117,100	IL6	Guinea pig
SEA079Hu	ELISA Kit for Interleukin 6 (IL6)	96Tests	114,000	79,800	IL6	Human
SEA079Mu	ELISA Kit for Interleukin 6 (IL6)	96Tests	117,200	82,000	IL6	Mouse
SEA079Ov	ELISA Kit for Interleukin 6 (IL6)	96Tests	188,400	131,800	IL6	Ovine
SEA079Po	ELISA Kit for Interleukin 6 (IL6)	96Tests	134,400	94,000	IL6	Porcine
SEA079Ra	ELISA Kit for Interleukin 6 (IL6)	96Tests	104,200	72,900	IL6	Rat
SEA079Rb	ELISA Kit for Interleukin 6 (IL6)	96Tests	159,100	111,300	IL6	Rabbit
SEA079Si	ELISA Kit for Interleukin 6 (IL6)	96Tests	188,400	131,800	IL6	Rhesus monkey
SEA087Ca	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	104,200	72,900	MCP-1	Canine
SEA087Hu	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	146,500	102,500	MCP-1	Human
SEA087Mi	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	150,700	105,400	MCP-1	データシートご確認下さい。
SEA087Mu	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	159,100	111,300	MCP-1	Mouse

Cloud-Clone社「新型コロナウイルス感染研究ターゲットELISA KIT」30%FFキャンペーン

対象製品価格表（2020年10/1～12/25）



Cloud-Clone Corp.

製品コード	品名	容量	希望納入価格 (円)	キャンペーン価格 (円)	ターゲット 因子	測定動物種
SEA087Po	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	175,800	123,000	MCP-1	Porcine
SEA087Ra	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	159,100	111,300	MCP-1	Rat
SEA087Rb	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	159,100	111,300	MCP-1	Rabbit
SEA087Si	ELISA Kit for Monocyte Chemotactic Protein 1 (MCP1)	96Tests	188,400	131,800	MCP-1	Rhesus monkey
SEA092Ca	ELISA Kit for Macrophage Inflammatory Protein 1 Alpha (MIP1a)	96Tests	167,400	117,100	MIP-1A	Canine
SEA092Hu	ELISA Kit for Macrophage Inflammatory Protein 1 Alpha (MIP1a)	96Tests	108,900	76,200	MIP-1A	Human
SEA092Mu	ELISA Kit for Macrophage Inflammatory Protein 1 Alpha (MIP1a)	96Tests	92,600	64,800	MIP-1A	Mouse
SEA092Ra	ELISA Kit for Macrophage Inflammatory Protein 1 Alpha (MIP1a)	96Tests	123,700	86,500	MIP-1A	Rat
SEA116Gu	ELISA Kit for Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES)	96Tests	167,400	117,100	CCL5	Guinea pig
SEA116Hu	ELISA Kit for Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES)	96Tests	108,900	76,200	CCL5	Human
SEA116Mu	ELISA Kit for Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES)	96Tests	92,600	64,800	CCL5	Mouse
SEA116Po	ELISA Kit for Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES)	96Tests	175,800	123,000	CCL5	Porcine
SEA116Ra	ELISA Kit for Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES)	96Tests	108,900	76,200	CCL5	Rat
SEA133Bo	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	156,300	109,400	TNF α	Bovine
SEA133Ca	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	130,200	91,100	TNF α	Canine
SEA133Cp	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	180,000	126,000	TNF α	Caprine
SEA133Eq	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	180,000	126,000	TNF α	Equine
SEA133Ga	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	159,100	111,300	TNF α	Chicken
SEA133Gu	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	130,200	91,100	TNF α	Guinea pig
SEA133Hu	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	104,200	72,900	TNF α	Human
SEA133Mi	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	111,600	78,100	TNF α	データシートご確認下さい。
SEA133Mu	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	117,200	82,000	TNF α	Mouse
SEA133Ov	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	188,400	131,800	TNF α	Ovine
SEA133Po	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	156,300	109,400	TNF α	Porcine
SEA133Ra	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	123,700	86,500	TNF α	Rat
SEA133Rb	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	123,700	86,500	TNF α	Rabbit
SEA133Si	ELISA Kit for Tumor Necrosis Factor Alpha (TNFa)	96Tests	104,200	72,900	TNF α	Rhesus monkey
SEA222Bo	ELISA Kit for Interferon Beta (IFNb)	96Tests	175,800	123,000	IFN β	Bovine
SEA222Ca	ELISA Kit for Interferon Beta (IFNb)	96Tests	167,400	117,100	IFN β	Canine
SEA222Eq	ELISA Kit for Interferon Beta (IFNb)	96Tests	180,000	126,000	IFN β	Equine
SEA222Ga	ELISA Kit for Interferon Beta (IFNb)	96Tests	159,100	111,300	IFN β	Chicken
SEA222Hu	ELISA Kit for Interferon Beta (IFNb)	96Tests	130,200	91,100	IFN β	Human
SEA222Mu	ELISA Kit for Interferon Beta (IFNb)	96Tests	150,700	105,400	IFN β	Mouse
SEA222Po	ELISA Kit for Interferon Beta (IFNb)	96Tests	175,800	123,000	IFN β	Porcine
SEA222Ra	ELISA Kit for Interferon Beta (IFNb)	96Tests	159,100	111,300	IFN β	Rat
SEA222Rb	ELISA Kit for Interferon Beta (IFNb)	96Tests	159,100	111,300	IFN β	Rabbit
SEA222Si	ELISA Kit for Interferon Beta (IFNb)	96Tests	188,400	131,800	IFN β	Rhesus monkey
SEA371Bo	ELISA Kit for Interferon Gamma Induced Protein 10kDa (IP10)	96Tests	175,800	123,000	CXCL10	Bovine
SEA371Ca	ELISA Kit for Interferon Gamma Induced Protein 10kDa (IP10)	96Tests	167,400	117,100	CXCL10	Canine
SEA371Hu	ELISA Kit for Interferon Gamma Induced Protein 10kDa (IP10)	96Tests	115,800	81,000	CXCL10	Human
SEA371Mu	ELISA Kit for Interferon Gamma Induced Protein 10kDa (IP10)	96Tests	102,100	71,400	CXCL10	Mouse
SEA371Ra	ELISA Kit for Interferon Gamma Induced Protein 10kDa (IP10)	96Tests	159,100	111,300	CXCL10	Rat
SEA821Bo	ELISA Kit for C Reactive Protein (CRP)	96Tests	175,800	123,000	CRP	Bovine
SEA821Ga	ELISA Kit for C Reactive Protein (CRP)	96Tests	159,100	111,300	CRP	Chicken
SEA821Gu	ELISA Kit for C Reactive Protein (CRP)	96Tests	167,400	117,100	CRP	Guinea pig
SEA821Hu	ELISA Kit for C Reactive Protein (CRP)	96Tests	111,100	77,700	CRP	Human
SEA821Mu	ELISA Kit for C Reactive Protein (CRP)	96Tests	115,800	81,000	CRP	Mouse
SEA821Po	ELISA Kit for C Reactive Protein (CRP)	96Tests	175,800	123,000	CRP	Porcine
SEA821Ra	ELISA Kit for C Reactive Protein (CRP)	96Tests	123,700	86,500	CRP	Rat
SEA821Rb	ELISA Kit for C Reactive Protein (CRP)	96Tests	159,100	111,300	CRP	Rabbit
SEA859Hu	ELISA Kit for Mannose Associated Serine Protease 2 (MASP2)	96Tests	146,500	102,500	MASP-2	Human

Cloud-Clone社「新型コロナウイルス感染研究ターゲットELISA KIT」30%FFキャンペーン
対象製品価格表（2020年10/1～12/25）



製品コード	品名	容量	希望納入価格 (円)	キャンペーン価格 (円)	ターゲット 因子	測定動物種
SEB589Hu	ELISA Kit for Interferon Regulatory Factor 3 (IRF3)	96Tests	154,700	108,200	IRF3	Human
SEB589Mu	ELISA Kit for Interferon Regulatory Factor 3 (IRF3)	96Tests	159,100	111,300	IRF3	Mouse
SEB740Hu	ELISA Kit for Signal Transducer And Activator Of Transcription 1 (STAT1)	96Tests	162,800	113,900	STAT1	Human
SEB740Mi	ELISA Kit for Signal Transducer And Activator Of Transcription 1 (STAT1)	96Tests	159,100	111,300	STAT1	データシートご確認下さい。
SEB743Hu	ELISA Kit for Signal Transducer And Activator Of Transcription 3 (STAT3)	96Tests	162,800	113,900	STAT3	Human
SEB743Mi	ELISA Kit for Signal Transducer And Activator Of Transcription 3 (STAT3)	96Tests	159,100	111,300	STAT3	データシートご確認下さい。
SEB796Hu	ELISA Kit for Signal Transducer And Activator Of Transcription 2 (STAT2)	96Tests	154,700	108,200	STAT2	Human
SEB796Ra	ELISA Kit for Signal Transducer And Activator Of Transcription 2 (STAT2)	96Tests	167,900	117,500	STAT2	Rat
SEB824Ga	ELISA Kit for Nuclear Factor Kappa B (NFkB)	96Tests	167,900	117,500	NF-κB	Chicken
SEB824Hu	ELISA Kit for Nuclear Factor Kappa B (NFkB)	96Tests	154,700	108,200	NF-κB	Human
SEB824Mu	ELISA Kit for Nuclear Factor Kappa B (NFkB)	96Tests	134,000	93,800	NF-κB	Mouse
SEB824Ra	ELISA Kit for Nuclear Factor Kappa B (NFkB)	96Tests	167,900	117,500	NF-κB	Rat
SEB950Hu	ELISA Kit for Toll Like Receptor 7 (TLR7)	96Tests	154,700	108,200	TLR7	Human
SEB950Mu	ELISA Kit for Toll Like Receptor 7 (TLR7)	96Tests	159,100	111,300	TLR7	Mouse
SEB989Hu	ELISA Kit for Toll Like Receptor 3 (TLR3)	96Tests	154,700	108,200	TLR3	Human
SEB989Mu	ELISA Kit for Toll Like Receptor 3 (TLR3)	96Tests	159,100	111,300	TLR3	Mouse
SEC795Hu	ELISA Kit for Transmembrane Protease, Serine 2 (TMPRSS2)	96Tests	162,800	113,900	TMPRSS2	Human
SEC795Mu	ELISA Kit for Transmembrane Protease, Serine 2 (TMPRSS2)	96Tests	167,400	117,100	TMPRSS2	Mouse
SEH061Hu	ELISA Kit for TANK Binding Kinase 1 (TBK1)	96Tests	162,800	113,900	TBK1	Human

- 掲載されている試薬は、試験・研究の目的にのみ使用されるものであり、「医薬品」、「食品」、「家庭用品」などとしては使用できません。
- 希望納入価格には消費税等が含まれておりません。
- 本キャンペーンは富士フイルム和光純薬株式会社の代理店・特約店から購入した場合のみ適用されます。
- 最新の保存条件・法規情報は当社、製品検索サイト(<https://labchem-wako.fujifilm.com/jp/index.html>)をご確認ください。