

富士フイルム和光純薬 抗体お試しキャンペーン

当社は汎用試薬だけでなく、神経科学領域やがん領域を中心に様々な抗体をラインアップしております。特にミクログリアマーカーの抗Iba1抗体は世界中の研究者に使用されています。

「市販の抗体は高いのに使えるか分からない…」
そんな研究者の悩みに応えるべく、抗体お試しキャンペーンを実施します。
この機会に是非、富士フイルム和光純薬の抗体をお試しください。

期間

2021年2月1日(月)～3月31日(水)

キャンペーン概要

① パンフレット掲載抗体が**最大50%OFF**

期間中、パンフレット掲載の抗体を最大50%OFFでご提供します。

② 小包装抗体を最大3種類 **無償提供** (初回購入者限定)


サンプルコードが記載されている抗体は10 μ L包装を無償提供します。
申込方法の詳細はパンフレット裏面をご覧ください。

注目の抗体

ロングセラーの抗Iba1抗体に橙色標識が登場

Anti Iba1, Rabbit, SPICA Dye™ 568-conjugated

ミクログリアマーカーIba1の抗体を新規耐光性蛍光色素SPICA Dye™ 568で標識。二次抗体処理を省略することができます。




ラット
ラット大脳皮質凍結切片の染色像

ヤギポリクローナル 抗Iba1抗体

Anti Iba1, Goat

ヤギポリクローナル 抗Iba1抗体はマウス、ラットサンプルの免疫組織染色およびウエスタンブロッティングに適用可能です。



ラット
ラット大脳皮質凍結切片の染色像

うつ病マーカー候補の注目分子BDNFを検出

Anti BDNF, Monoclonal Antibody

脳由来神経栄養因子 (BDNF) は神経の発達、維持に重要な分子です。

特にうつ病においてはKaregeらがうつ病患者の血清中におけるBDNFが健常者に比べて低値であることを初めて報告して以降¹⁾、バイオマーカーとしての利用が期待されています。

1) Karege et al., Decreased Serum Brain-Derived Neurotrophic Factor Levels in Major Depressed Patients. *Psychiatry Res* 109 (2) : 143-148 (2002).

キャンペーン対象製品

標的抗原	コードNo.	製品名 [クローンNo.]	容量	免疫動物	交差性	適用	希望納入価格 (円)	キャンペーン価格 (円)	サンプルコード
グリア細胞									
FGF9	014-27383	Anti Human FGF9, Monoclonal Antibody [GAF150-59]	50μL	Mouse	Human	ELISA, Neu	35,000	17,500	A
GFAP	018-27283	Anti GFAP, Monoclonal Antibody [MO389]	50μL	Mouse	Human	ICC, WB	30,000	21,000	-
Phosphorylated GFAP S8	018-27303	Anti Human Phosphorylated GFAP S8, Monoclonal Antibody [YC10]	50μL	Mouse	Human	ICC, WB	30,000	15,000	B
Phosphorylated GFAP S13	012-27323	Anti Human Phosphorylated GFAP S13, Monoclonal Antibody [KT13]	50μL	Mouse	Human	ICC, WB	30,000	15,000	C
Iba1	016-26721	Anti Iba1, Monoclonal Antibody [NCNP24]	50μL	Mouse	Marmoset, Mouse, Rat	IHC(F)	45,000	31,500	D
Iba1	013-27593	Anti Iba1, Monoclonal Antibody [NCNP27]	50μL	Mouse	Human	IHC(P)	40,000	28,000	-
Iba1	011-27991	Anti Iba1, Goat [Polyclonal]	100μL	Goat	Mouse, Rat	IHC(F), IHC(P), WB	45,000	31,500	-
Iba1	015-28011	Anti Iba1, Rabbit, SPICA Dye™ 568-conjugated [Polyclonal]	100μL	Rabbit	Mouse, Rat	IHC(F)	50,000	25,000	-
アルツハイマー病									
Amyloid Precursor Protein	010-27243	Anti Amyloid Precursor Protein, Monoclonal Antibody [3E9]	50μL	Mouse	Human, Mouse	WB	30,000	15,000	E
Amyloid β	010-26883	Anti Amyloid β, Monoclonal Antibody [BNT77]	50μL	Mouse	Human, Mouse, Rat	ELISA, IHC, IP,	35,000	24,500	F
Amyloid β	012-27781	Anti Amyloid β, Monoclonal Antibody [2C8]	50μL	Mouse	Human, Mouse	IHC, WB	35,000	24,500	-
Amyloid β	013-26873	Anti Human Amyloid β, Monoclonal Antibody [BAN50]	50μL	Mouse	Human	ELISA, ICC, IHC, IP, WB	35,000	24,500	-
Amyloid β40	014-26923	Anti Amyloid β40, Monoclonal Antibody [BA27]	50μL	Mouse	Human, Mouse, Rat	ELISA, IHC, WB	35,000	24,500	-
Amyloid β42 (43)	010-26903	Anti Amyloid β42(43), Monoclonal Antibody [BC05]	50μL	Mouse	Human, Mouse, Rat	ELISA, IHC, WB	35,000	24,500	-
Tau	017-26893	Anti Human/Mouse/Rat Tau, Rat Monoclonal Antibody [RTM38]	50μL	Rat	Human, Mouse, Rat	ICC, IHC, WB	35,000	24,500	-
Tau (Human specific)	015-26953	Anti Human Tau, Rat Monoclonal Antibody [RTM49]	50μL	Rat	Human	ICC, IHC, WB	35,000	24,500	G
Tau (Mouse specific)	012-26963	Anti Mouse Tau, Rat Monoclonal Antibody [RTM47]	50μL	Rat	Mouse	ICC, IHC, WB	35,000	24,500	H
2N-Tau	013-27353	Anti 2N-Tau, Rat Monoclonal Antibody [2C2]	50μL	Rat	Human	IHC, WB	30,000	15,000	I
3R-Tau	016-26581	Anti 3R-Tau, Rat Monoclonal Antibody [2A1-1F4]	50μL	Rat	Human	IP, WB	30,000	21,000	J
4R-Tau	013-26591	Anti 4R-Tau, Monoclonal Antibody [3E8-1A6]	50μL	Mouse	Human	IP, WB	30,000	21,000	K
Phosphorylated Tau S422	012-27683	Anti Phosphorylated Tau S422, Monoclonal Antibody [AP422]	50μL	Mouse	Human, Mouse	IHC, WB	50,000	35,000	L
Phosphorylated Tau T181	016-26601	Anti Phosphorylated Tau T181, Rat Monoclonal Antibody [2E2-A6]	50μL	Rat	Human	IHC, WB	30,000	21,000	M
Apolipoprotein E4	014-27263	Anti Human Apolipoprotein E4, Monoclonal Antibody [1F9]	50μL	Mouse	Human	IHC, WB	30,000	21,000	N
CRMP2	014-24821	Anti CRMP2, Monoclonal Antibody [9F]	100μL	Mouse	Chicken, Human, Mouse, Rat	ELISA, ICC, IHC, WB	40,000	28,000	-

【アプリケーション略号一覧】

ELISA：酵素結合免疫吸着測定法 / FC：フローサイトメトリー / ICC：免疫細胞染色 / IHC：免疫組織染色 / IHC(F)：免疫組織染色(凍結切片) / IHC(P)：免疫組織染色(パラフィン切片) / IP：免疫沈降 / WB：ウエスタンブロッティング / Neu：中和活性あり / ADCC：抗体依存性細胞傷害あり / CDC：補体依存性細胞傷害あり

標的抗原	コードNo.	製品名 [クローンNo.]	容量	免疫動物	交差性	適用	希望納入価格 (円)	キャンペーン価格 (円)	サンプルコード
パーキンソン病									
Parkin	016-27343	Anti Parkin, Monoclonal Antibody [Par6]	50μL	Mouse	Human, Mouse, Rat	WB	30,000	15,000	-
Synaptophysin	019-27333	Anti Synaptophysin, Monoclonal Antibody [171B5]	50μL	Mouse	Human, Mouse, Rat	IHC, WB	30,000	21,000	-
Phosphorylated α-Synuclein	010-26481	Anti Phosphorylated α-Synuclein, Monoclonal Antibody, Biotin-conjugated [pSyn#64]	100μL	Mouse	Human, Mouse, Rat	ICC, IHC	50,000	35,000	-
神経栄養因子									
BDNF	017-28071	Anti BDNF, Monoclonal Antibody [2D7]	100μL	Mouse	Human, Mouse	ELISA, IP	50,000	25,000	-
BDNF	014-28081	Anti BDNF, Monoclonal Antibody [3C8]	100μL	Mouse	Human, Mouse	ELISA, IP	50,000	25,000	-
NGF	014-27143	Anti NGF, Monoclonal Antibody [NGFA-133]	50μL	Mouse	Human, Mouse, Rat	ELISA, Neu	35,000	17,500	-
NT-3	011-27153	Anti NT-3, Monoclonal Antibody [3W3]	50μL	Mouse	Human, Mouse, Rat	ELISA, IHC(F), Neu	35,000	17,500	-
神経発生・成長									
Drebrin	011-27273	Anti Drebrin, Monoclonal Antibody [M2F6]	50μL	Mouse	Human, Mouse, Rat	ICC, IHC, WB	30,000	21,000	-
Phosphorylated GAP-43 S96	017-25391	Anti Phosphorylated GAP-43 S96, Monoclonal Antibody [16-4C]	100μL	Mouse	Mouse, Rat	IHC, WB	48,500	33,950	-
Phosphorylated GAP-43 S96	010-25401	Anti Phosphorylated GAP-43 S96, Monoclonal Antibody [18-10H-9H]	100μL	Mouse	Mouse, Rat	IHC, WB	48,500	33,950	-
Phosphorylated GAP-43 T172	017-25411	Anti Phosphorylated GAP-43 T172, Monoclonal Antibody [19-9A]	100μL	Mouse	Human, Mouse, Rat	IHC	48,500	33,950	-
Pax6	015-27293	Anti Pax6, Rabbit [Polyclonal]	50μL	Rabbit	Chicken, Mouse, Rat	IHC, WB	30,000	21,000	-
Nestin	012-26843	Anti Mouse Nestin, Rat Monoclonal Antibody [7A3]	50μL	Rat	Mouse	ICC, IHC	30,000	21,000	-
血液脳関門									
Claudin-5	014-28101	Anti Claudin-5, Rat Monoclonal Antibody [R9]	500μL	Rat	Human, Monkey	FC, Neu	70,000	49,000	-
Claudin-5	011-28091	Anti Claudin-5, Monoclonal Antibody [M48]	500μL	Mouse	Human, Monkey	FC, Neu	70,000	49,000	-
神経腫									
IDH1-R132H mutation	018-24081	Anti IDH1-R132H, Monoclonal Antibody [HMab-1]	100μg	Mouse	Human	ELISA, IHC, WB	36,400	25,480	-
IDH1-R132H mutation	013-26851	Anti IDH1-R132H, Monoclonal Antibody [HMab-2]	100μg	Mouse	Human	ELISA, IHC, WB	40,000	28,000	-
Mutated IDH1/2	015-25691	Anti Mutated IDH1/2, Monoclonal Antibody [MsMab-1]	100μg	Mouse	Human	ELISA, IHC, WB	36,400	25,480	-
TERT	010-26861	Anti TERT, Monoclonal Antibody [TMab-6]	100μg	Mouse	Human	ELISA, IHC, WB	40,000	28,000	-
がん									
HB-EGF	019-27193	Anti HB-EGF, Monoclonal Antibody [Y-073]	50μL	Mouse	Human, Mouse, Rat	ELISA, Neu	35,000	17,500	-
Podoplanin	018-24101	Anti Human Podoplanin, Monoclonal Antibody [NZ-1.2]	100μg	Rat	Human	ELISA, FC, IHC, IP, WB	36,400	25,480	-
Podoplanin	013-27071	Anti Human Podoplanin, Monoclonal Antibody [LpMab-10]	100μg	Mouse	Human	ELISA, FC, ICC, IHC, WB	34,000	23,800	-
Podoplanin	014-27001	Anti Human Podoplanin, Monoclonal Antibody [LpMab-12]	100μg	Mouse	Human	ELISA, FC, IHC, Neu, WB	34,000	23,800	-
Podoplanin	011-27011	Anti Human Podoplanin, Monoclonal Antibody [LpMab-13]	100μg	Mouse	Human	ELISA, FC, IHC, WB	34,000	23,800	-
Podoplanin	018-27021	Anti Human Podoplanin, Monoclonal Antibody [LpMab-17]	100μg	Mouse	Human	ELISA, FC, IHC, WB	34,000	23,800	-
Podoplanin	017-27471	Anti Human Podoplanin, Monoclonal Antibody [LpMab-19]	100μg	Mouse	Human	ELISA, FC, IHC	34,000	23,800	-
Podoplanin	014-27481	Anti Human Podoplanin, Monoclonal Antibody [LpMab-21]	100μg	Mouse	Human	ADCC, CDC, ELISA, FC, IHC	34,000	23,800	-
Podoplanin	012-27041	Anti Rabbit Podoplanin, Monoclonal Antibody [PMab-32]	100μg	Mouse	Rabbit	ELISA, FC, IHC, WB	34,000	23,800	-
Podoplanin	015-24111	Anti Mouse Podoplanin, Monoclonal Antibody [PMab-1]	100μg	Rat	Mouse	ELISA, FC, IHC, IP, WB	36,400	25,480	-
Podoplanin	017-27091	Anti Dog Podoplanin, Monoclonal Antibody [PMab-38]	100μg	Mouse	Dog	ELISA, FC, IHC, WB	34,000	23,800	-
Podocalyxin	011-27491	Anti Human Podocalyxin, Monoclonal Antibody [PcMab-47b]	100μg	Mouse	Human	ELISA, FC, IHC	34,000	23,800	-
線維化									
CTGF Module2	015-27433	Anti CTGF Module2, Monoclonal Antibody [2-3]	50μL	Mouse	Human, Mouse	Neu	45,000	31,500	O
CTGF Module3	012-27443	Anti CTGF Module3, Monoclonal Antibody [3-54]	50μL	Mouse	Human, Mouse	Neu	45,000	31,500	P
CTGF Module4	019-27453	Anti CTGF Module4, Monoclonal Antibody [4-69]	50μL	Mouse	Human, Mouse	Neu	45,000	31,500	Q
Integrin α8	019-28053	Anti Integrin α8, Monoclonal Antibody, Chicken / Mouse IgG1-κ Chimera [YZ26]	50μL	Chicken Mouse※	Human, Mouse, Rat	FC, ICC, Neu	35,000	17,500	-
Integrin α9	016-28063	Anti Integrin α9, Monoclonal Antibody, Chicken / Mouse IgG1-κ Chimera [YA4]	50μL	Chicken Mouse※	Human, Mouse, Rat	FC, ICC, Neu	35,000	17,500	-

※ 抗Integrin α8/α9 抗体は可変領域：ニフトリ、定常領域：マウスIgG1のキメラ抗体です。二次抗体には抗マウス抗体をおすすめいたします。

抗体お試しキャンペーン 無償提供申込用紙

Step1: ご希望の無償提供製品に✓をつけてください。(最大3種類まで選択可)

- 抗体の無償提供は、対象抗体を始めて使用する方に限ります。
- 1研究室・1グループに1回のみのご提供とさせていただきます。また1製品につき1本、最大3種類のご提供となります。
- 定員に達した場合、申込みされても無償提供できない可能性がございます。(最新の配布状況は当社Webサイトにてご確認ください)

サンプルコード	製品名 [クローンNo.]	申込	サンプルコード	製品名 [クローンNo.]	申込
A	Anti Human FGF9, Monoclonal Antibody [GAF150-59]		I	Anti 2N-Tau, Rat Monoclonal Antibody [2C2]	
B	Anti Human Phosphorylated GFAP S8, Monoclonal Antibody [YC10]		J	Anti 3R-Tau, Rat Monoclonal Antibody [2A1-1F4]	
C	Anti Human Phosphorylated GFAP S13, Monoclonal Antibody [KT13]		K	Anti 4R-Tau, Monoclonal Antibody [3E8-1A6]	
D	Anti Iba1, Monoclonal Antibody [NCNP24]		L	Anti Phosphorylated Tau S422, Monoclonal Antibody [AP422]	
E	Anti Amyloid Precursor Protein, Monoclonal Antibody [3E9]		M	Anti Phosphorylated Tau T181, Rat Monoclonal Antibody [2E2-A6]	
F	Anti Amyloid β , Monoclonal Antibody [BNT77]		N	Anti Human Apolipoprotein E4, Monoclonal Antibody [1F9]	
G	Anti Human Tau, Rat Monoclonal Antibody [RTM49]		O	Anti CTGF Module2, Monoclonal Antibody [2-3]	
H	Anti Mouse Tau, Rat Monoclonal Antibody [RTM47]		P	Anti CTGF Module3, Monoclonal Antibody [3-54]	
			Q	Anti CTGF Module4, Monoclonal Antibody [4-69]	

Step2: 各種キャンペーンや新製品をはじめ最新情報満載のメールマガジン配信を希望されますか。

- ☐ 配信を希望する。 ☐ 配信を希望しない。 ☐ 既に登録している。

Step3: お客様の情報をご記入ください。

所属			
氏名		販売代理店	
E-mail		TEL	

Step4: 本申込用紙をFAX(03-3270-8547)でお送りいただくか、当社営業もしくは販売代理店にお渡しください。

- 本文に収載しております試薬は、試験・研究の目的にのみ使用されるもので、「医薬品」、「食品」、「生活用品」などとして使用できません。
- 希望納入価格、キャンペーン価格には消費税等が含まれておりません。
- 本キャンペーンは富士フイルム和光純薬(株)の代理店・特約店から購入した場合のみ適用されます。

富士フイルム 和光純薬株式会社

本 社 〒540-8605 大阪市中央区道修町三丁目1番2号 TEL 06-6203-3741 (代表)
東京本店 〒103-0023 東京都中央区日本橋本町二丁目4番1号 TEL 03-3270-8571 (代表)

- 九州営業所 ● 中国営業所
● 東海営業所 ● 横浜営業所
● 筑波営業所 ● 東北営業所
● 北海道営業所
-  フリーダイヤル 0120-052-099
試薬URL: <https://labchem-wako.fujifilm.com>

■ FUJIFILM Wako Chemicals U.S.A. Corporation
1600 Bellwood Road, Richmond, VA 23237, USA
TEL: +1-804-714-1920 FAX: +1-804-271-7791

■ 富士フイルム和光(香港)有限公司
Room 1111, 11/F, International Trade Centre,
11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong
TEL: +852-2799-9019 FAX: +852-2799-9808

■ FUJIFILM Wako Chemicals Europe GmbH
Fuggerstr 12, 41468 Neuss, Germany
TEL: +49-2131-311-0 FAX: +49-2131-311-100

■ 富士フイルム和光(広州)貿易有限公司
广州市越秀区先烈中路69号东山广场30楼3002-3003室
TEL: +86-20-8732-6381 (广州)
TEL: +86-21-6288-4751 (上海)
TEL: +86-10-6413-6388 (北京)