

●UnyLinker™ Supports

担体	リンカー	ローディング量 ( $\mu\text{mol/g}$ )	品名	メーカーコード	容量	コードNo.	希望納入価格(円)	キャンペーン価格(円)
NitroPhase™	UnyLinker™	50	NitroPhase(TM) UnyLinker(TM)50	10-00-03-050	1g	553-31631	26,200	20,960
					5g	559-31633	112,100	89,680
					10g	557-31634	186,200	148,960
					50g	553-31636	931,000	744,800
					100g	551-31637	1,482,000	1,185,600
		80	NitroPhase(TM) UnyLinker(TM)80	10-00-03-080	1g	550-31641	26,200	20,960
					5g	556-31643	112,100	89,680
					10g	554-31644	186,200	148,960
					50g	550-31646	931,000	744,800
					100g	558-31647	1,482,000	1,185,600
		100	NitroPhase(TM) UnyLinker(TM)100	10-00-03-100	1g	557-31651	26,200	20,960
					5g	553-31653	112,100	89,680
					10g	551-31654	186,200	148,960
					50g	557-31656	931,000	744,800
					100g	555-31657	1,482,000	1,185,600
		200	NitroPhase(TM) UnyLinker(TM)200	10-00-03-200	1g	554-31661	26,200	20,960
					5g	550-31663	112,100	89,680
					10g	558-31664	148,200	118,560
					50g	554-31666	741,000	592,800
					100g	552-31667	1,102,000	881,600
NitroPhase™ HL	UnyLinker™	150	NitroPhase(TM) HL UnyLinker(TM)150	20-00-03-150	1g	551-31671	26,200	20,960
					5g	557-31673	112,100	89,680
					10g	555-31674	186,200	148,960
					50g	551-31676	931,000	744,800
					100g	559-31677	1,482,000	1,185,600
		200	NitroPhase(TM) HL UnyLinker(TM)200	20-00-03-200	1g	558-31681	26,200	20,960
					5g	554-31683	112,100	89,680
					10g	552-31684	186,200	148,960
					50g	558-31686	931,000	744,800
					100g	556-31687	1,482,000	1,185,600
		250	NitroPhase(TM) HL UnyLinker(TM)250	20-00-03-250	1g	555-31691	26,200	20,960
					5g	551-31693	112,100	89,680
					10g	559-31694	186,200	148,960
					50g	555-31696	931,000	744,800
					100g	553-31697	1,482,000	1,185,600
		300	NitroPhase(TM) HL UnyLinker(TM)300	20-00-03-300	1g	558-31701	26,200	20,960
					5g	554-31703	112,100	89,680
					10g	552-31704	186,200	148,960
					50g	558-31706	931,000	744,800
					100g	556-31707	1,482,000	1,185,600
		350	NitroPhase(TM) HL UnyLinker(TM)350	20-00-03-350	1g	555-31711	26,200	20,960
					5g	551-31713	112,100	89,680
					10g	559-31714	148,200	118,560
					50g	555-31716	741,000	592,800
					100g	553-31717	1,102,000	881,600
		400	NitroPhase(TM) HL UnyLinker(TM)400	20-00-03-400	1g	552-31721	26,200	20,960
					5g	558-31723	112,100	89,680
					10g	556-31724	186,200	148,960
					50g	552-31726	931,000	744,800
					100g	550-31727	1,482,000	1,185,600

●2'-OMe Supports

担体	リンカー	ローディング量 ( $\mu\text{mol/g}$ )	品名	メーカーコード	容量	コードNo.	希望納入価格(円)	キャンペーン価格(円)
NitroPhase™ HL	2'-OMe A(Bz)	200	NitroPhase(TM) HL 2'OMeA(bz)200	24-12-00-200	1g	557-32511	36,100	28,880
					5g	553-32513	150,100	120,080
					10g	551-32514	224,200	179,360
					50g	559-32515	1,121,000	896,800
					100g	557-32516	1,862,000	1,489,600
		250	NitroPhase(TM) HL 2'OMeA(bz)250	24-12-00-250	1g	554-32521	36,100	28,880
					5g	550-32523	150,100	120,080
					10g	558-32524	224,200	179,360
					50g	556-32525	1,121,000	896,800
					100g	554-32526	1,862,000	1,489,600
		300	NitroPhase(TM) HL 2'OMeA(bz)300	24-12-00-300	1g	557-54731	36,100	28,880
					5g	553-54733	150,100	120,080
					10g	551-54734	224,200	179,360
					50g	557-54736	1,121,000	896,800
					100g	555-54732	1,862,000	1,489,600
		350	NitroPhase(TM) HL 2'OMeA(bz)350	24-12-00-350	1g	551-32531	36,100	28,880
					5g	557-32533	150,100	120,080
					10g	555-32534	224,200	179,360
					50g	553-32535	1,121,000	896,800
					100g	551-32536	1,862,000	1,489,600
NitroPhase™ HL	2'-OMe U	200	NitroPhase(TM) HL 2'OMeU200	24-50-00-200	1g	553-32611	36,100	28,880
					5g	559-32613	150,100	120,080
					10g	557-32614	224,200	179,360
					50g	555-32615	1,121,000	896,800
					100g	553-32616	1,862,000	1,489,600
		250	NitroPhase(TM) HL 2'OMeU250	24-50-00-250	1g	550-32621	36,100	28,880
					5g	556-32623	150,100	120,080
					10g	554-32624	224,200	179,360
					50g	552-32625	1,121,000	896,800
					100g	550-32626	1,862,000	1,489,600
		300	NitroPhase(TM) HL 2'OMeU300	24-50-00-300	1g	558-54761	36,100	28,880
					5g	554-54763	150,100	120,080
					10g	552-54764	224,200	179,360
					50g	558-54766	1,121,000	896,800
					100g	556-54762	1,862,000	1,489,600
		350	NitroPhase(TM) HL 2'OMeU350	24-50-00-350	1g	557-32631	36,100	28,880
					5g	553-32633	150,100	120,080
					10g	551-32634	224,200	179,360
					50g	559-32635	1,121,000	896,800
					100g	557-32636	1,862,000	1,489,600

NittoPhase™ HL	2'-OMe G(iBu)	200	NittoPhase(TM) HL 2'OMeG(ibu)200	24-43-00-200	1g	556-32581	36,100	28,880
					5g	552-32583	150,100	120,080
					10g	550-32584	224,200	179,360
					50g	558-32585	1,121,000	896,800
					100g	556-32586	1,862,000	1,489,600
		250	NittoPhase(TM) HL 2'OMeG(ibu)250	24-43-00-250	1g	553-32591	36,100	28,880
					5g	559-32593	150,100	120,080
					10g	557-32594	224,200	179,360
					50g	555-32595	1,121,000	896,800
					100g	553-32596	1,862,000	1,489,600
		300	NittoPhase(TM) HL 2'OMeG(ibu)300	24-43-00-300	1g	551-54751	36,100	28,880
					5g	557-54753	150,100	120,080
					10g	555-54754	224,200	179,360
					50g	551-54756	1,121,000	896,800
					100g	559-54752	1,862,000	1,489,600
		350	NittoPhase(TM) HL 2'OMeG(ibu)350	24-43-00-350	1g	556-32601	36,100	28,880
					5g	552-32603	150,100	120,080
					10g	550-32604	224,200	179,360
					50g	558-32605	1,121,000	896,800
					100g	556-32606	1,862,000	1,489,600
NittoPhase™ HL	2'-OMe C(Ac)	200	NittoPhase(TM) HL 2'OMeC(ac)200	24-21-00-200	1g	558-32541	36,100	28,880
					5g	554-32543	150,100	120,080
					10g	552-32544	224,200	179,360
					50g	550-32545	1,121,000	896,800
					100g	558-32546	1,862,000	1,489,600
		250	NittoPhase(TM) HL 2'OMeC(ac)250	24-21-00-250	1g	555-32551	36,100	28,880
					5g	551-32553	150,100	120,080
					10g	559-32554	224,200	179,360
					50g	557-32555	1,121,000	896,800
					100g	555-32556	1,862,000	1,489,600
		300	NittoPhase(TM) HL 2'OMeC(ac)300	24-21-00-300	1g	554-54741	36,100	28,880
					5g	550-54743	150,100	120,080
					10g	558-54744	224,200	179,360
					50g	554-54746	1,121,000	896,800
					100g	552-54742	1,862,000	1,489,600
		350	NittoPhase(TM) HL 2'OMeC(ac)350	24-21-00-350	1g	559-32571	36,100	28,880
					5g	555-32573	150,100	120,080
					10g	553-32574	224,200	179,360
					50g	551-32575	1,121,000	896,800
					100g	559-32576	1,862,000	1,489,600

●DNA Supports

担体	リンカー	ローディング量 ( $\mu\text{mol/g}$ )	品名	メーカーコード	容量	コードNo.	希望納入価格(円)	キャンペーン価格(円)
NitroPhase™ HL	dA(Bz)	200	NitroPhase(TM) HL dA(bz)200	21-12-00-200	1g	557-54611	36,100	28,880
					5g	553-54613	150,100	120,080
					10g	551-54614	186,200	148,960
					50g	557-54616	931,000	744,800
					100g	555-54612	1,482,000	1,185,600
		250	NitroPhase(TM) HL dA(bz)250	21-12-00-250	1g	553-32351	36,100	28,880
					5g	559-32353	150,100	120,080
					10g	557-32354	186,200	148,960
					50g	555-32355	931,000	744,800
					100g	553-32356	1,482,000	1,185,600
		300	NitroPhase(TM) HL dA(bz)300	21-12-00-300	1g	554-54621	36,100	28,880
					5g	550-54623	150,100	120,080
					10g	558-54624	186,200	148,960
					50g	554-54626	931,000	744,800
					100g	552-54622	1,482,000	1,185,600
		350	NitroPhase(TM) HL dA(bz)350	21-12-00-350	1g	550-32361	36,100	28,880
					5g	556-32363	150,100	120,080
					10g	554-32364	186,200	148,960
					50g	552-32365	931,000	744,800
					100g	550-32366	1,482,000	1,185,600
NitroPhase™ HL	dT	200	NitroPhase(TM) HL T200	21-30-00-200	1g	555-54651	36,100	28,880
					5g	551-54653	150,100	120,080
					10g	559-54654	186,200	148,960
					50g	555-54656	931,000	744,800
					100g	553-54652	1,482,000	1,185,600
		250	NitroPhase(TM) HL T250	21-30-00-250	1g	551-32391	36,100	28,880
					5g	557-32393	150,100	120,080
					10g	555-32394	186,200	148,960
					50g	553-32395	931,000	744,800
					100g	551-32396	1,482,000	1,185,600
		300	NitroPhase(TM) HL T300	21-30-00-300	1g	552-54661	36,100	28,880
					5g	558-54663	150,100	120,080
					10g	556-54664	186,200	148,960
					50g	552-54666	931,000	744,800
					100g	550-54662	1,482,000	1,185,600
		350	NitroPhase(TM) HL T350	21-30-00-350	1g	554-32401	36,100	28,880
					5g	550-32403	150,100	120,080
					10g	558-32404	186,200	148,960
					50g	556-32405	931,000	744,800
					100g	554-32406	1,482,000	1,185,600

NittoPhase™ HL	dG(iBu)	200	NittoPhase(TM) HL dG(ibu)200	21-43-00-200	1g	559-54671	36,100	28,880
					5g	555-54673	150,100	120,080
					10g	553-54674	186,200	148,960
					50g	559-54676	931,000	744,800
					100g	557-54672	1,482,000	1,185,600
		250	NittoPhase(TM) HL dG(ibu)250	21-43-00-250	1g	551-32411	36,100	28,880
					5g	557-32413	150,100	120,080
					10g	555-32414	186,200	148,960
					50g	553-32415	931,000	744,800
					100g	551-32416	1,482,000	1,185,600
		300	NittoPhase(TM) HL dG(ibu)300	21-43-00-300	1g	556-54681	36,100	28,880
					5g	552-54683	150,100	120,080
					10g	550-54684	186,200	148,960
					50g	556-54686	931,000	744,800
					100g	554-54682	1,482,000	1,185,600
		350	NittoPhase(TM) HL dG(ibu)350	21-43-00-350	1g	558-32421	36,100	28,880
					5g	554-32423	150,100	120,080
					10g	552-32424	186,200	148,960
					50g	550-32425	931,000	744,800
					100g	558-32426	1,482,000	1,185,600
NittoPhase™ HL	dC(Ac)	200	NittoPhase(TM) HL dC(ac)200	21-21-00-200	1g	551-54631	36,100	28,880
					5g	557-54633	150,100	120,080
					10g	555-54634	186,200	148,960
					50g	551-54636	931,000	744,800
					100g	559-54632	1,482,000	1,185,600
		250	NittoPhase(TM) HL dC(ac)250	21-21-00-250	1g	557-32371	36,100	28,880
					5g	553-32373	150,100	120,080
					10g	551-32374	186,200	148,960
					50g	559-32375	931,000	744,800
					100g	557-32376	1,482,000	1,185,600
		300	NittoPhase(TM) HL dC(ac)300	21-21-00-300	1g	558-54641	36,100	28,880
					5g	554-54643	150,100	120,080
					10g	552-54644	186,200	148,960
					50g	558-54646	931,000	744,800
					100g	556-54642	1,482,000	1,185,600
		350	NittoPhase(TM) HL dC(ac)350	21-21-00-350	1g	554-32381	36,100	28,880
					5g	550-32383	150,100	120,080
					10g	558-32384	186,200	148,960
					50g	556-32385	931,000	744,800
					100g	554-32386	1,482,000	1,185,600

●RNA Supports

担体	リンカー	ローディング量 ( $\mu\text{mol/g}$ )	品名	メーカーコード	容量	コードNo.	希望納入価格(円)	キャンペーン価格(円)
NittoPhase™ HL	rA(Bz)	150	NittoPhase(TM) HL rA(bz)150	22-12-00-150	1g	555-32431	36,100	28,880
					5g	551-32433	152,000	121,600
					10g	559-32434	262,200	209,760
					50g	557-32435	1,311,000	1,048,800
					100g	555-32436	2,242,000	1,793,600
		200	NittoPhase(TM) HL rA(bz)200	22-12-00-200	1g	553-54691	36,100	28,880
					5g	559-54693	152,000	121,600
					10g	557-54694	262,200	209,760
					50g	553-54696	1,311,000	1,048,800
					100g	551-54692	2,242,000	1,793,600
		250	NittoPhase(TM) HL rA(bz)250	22-12-00-250	1g	552-32441	36,100	28,880
					5g	558-32443	152,000	121,600
					10g	556-32444	262,200	209,760
					50g	554-32445	1,311,000	1,048,800
					100g	552-32446	2,242,000	1,793,600
NittoPhase™ HL	rU	150	NittoPhase(TM) HL rU150	22-50-00-150	1g	557-32491	36,100	28,880
					5g	553-32493	152,000	121,600
					10g	551-32494	262,200	209,760
					50g	559-32495	1,311,000	1,048,800
					100g	557-32496	2,242,000	1,793,600
		200	NittoPhase(TM) HL rU200	22-50-00-200	1g	550-54721	36,100	28,880
					5g	556-54723	152,000	121,600
					10g	554-54724	262,200	209,760
					50g	550-54726	1,311,000	1,048,800
					100g	558-54722	2,242,000	1,793,600
		250	NittoPhase(TM) HL rU250	22-50-00-250	1g	550-32501	36,100	28,880
					5g	556-32503	152,000	121,600
					10g	554-32504	262,200	209,760
					50g	552-32505	1,311,000	1,048,800
					100g	550-32506	2,242,000	1,793,600
NittoPhase™ HL	rG(iBu)	150	NittoPhase(TM) HL rG(ibu)150	22-43-00-150	1g	553-32471	36,100	28,880
					5g	559-32473	152,000	121,600
					10g	557-32474	262,200	209,760
					50g	555-32475	1,311,000	1,048,800
					100g	553-32476	2,242,000	1,793,600
		200	NittoPhase(TM) HL rG(ibu)200	22-43-00-200	1g	553-54711	36,100	28,880
					5g	559-54713	152,000	121,600
					10g	557-54714	262,200	209,760
					50g	553-54716	1,311,000	1,048,800
					100g	551-54712	2,242,000	1,793,600
		250	NittoPhase(TM) HL rG(ibu)250	22-43-00-250	1g	550-32481	36,100	28,880
					5g	556-32483	152,000	121,600
					10g	554-32484	262,200	209,760
					50g	552-32485	1,311,000	1,048,800
					100g	550-32486	2,242,000	1,793,600
NittoPhase™ HL	rC(Ac)	150	NittoPhase(TM) HL rC(ac)150	22-21-00-150	1g	559-32451	36,100	28,880
					5g	555-32453	152,000	121,600
					10g	553-32454	262,200	209,760
					50g	551-32455	1,311,000	1,048,800
					100g	559-32456	2,242,000	1,793,600
		200	NittoPhase(TM) HL rC(ac)200	22-21-00-200	1g	556-54701	36,100	28,880
					5g	552-54703	152,000	121,600
					10g	550-54704	262,200	209,760
					50g	556-54706	1,311,000	1,048,800
					100g	554-54702	2,242,000	1,793,600
		250	NittoPhase(TM) HL rC(ac)250	22-21-00-250	1g	556-32461	36,100	28,880
					5g	552-32463	152,000	121,600
					10g	550-32464	262,200	209,760
					50g	558-32465	1,311,000	1,048,800
					100g	556-32466	2,242,000	1,793,600