

Code No.01-0170

SPECIFICATION

Aluminon

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|---------------------------------------|
| * Appearance | Dark reddish brown, crystals ~ powder |
| Solubility in water | to pass test |
| Loss on drying over silica gel in vacuum | max.5.0%(mass/mass) |
| Residue after ignition (as sulfate) | max.0.2%(mass/mass) |
| Sensitivity | to pass test |

JIS K 8011:2010

(*:Additional test performed by Wako but not required by JIS)