

Code No.01-0227

## SPECIFICATION

## 4-Aminoantipyrine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Pale yellow ~ pale yellowish brown, crystals ~
Assay (C <sub>11</sub> H <sub>13</sub> N <sub>3</sub> O)	crystalline powder min.99.0%(mass/mass)
Solubility in water	to pass test
Solubility in ethanol	to pass test
Melting point	107~110°C
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Suitability for phenol determination	to pass test

JIS K 8048:2021(\*:Additional test performed by Wako but not required by JIS)