

## SPECIFICATION

## Ammonia Solution

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (NH <sub>3</sub> )	25.0~27.9%(mass/mass)
Residue after evaporation	max.0.002%(mass/mass)
Carbonate (CO <sub>3</sub> )	max.0.002%(mass/mass)
Chloride (Cl)	max.0.3ppm(mass/mass)
Phosphate (PO <sub>4</sub> )	max.0.5ppm(mass/mass)
Silicate (as SiO <sub>2</sub> )	max.0.001%(mass/mass)
Sulfate (SO <sub>4</sub> )	max.2ppm(mass/mass)
Sulfide (S)	max.0.2ppm(mass/mass)
Sodium (Na)	max.1ppm(mass/mass)
Potassium (K)	max.1ppm(mass/mass)
Copper (Cu)	max.0.1ppm(mass/mass)
Magnesium (Mg)	max.1ppm(mass/mass)
Calcium (Ca)	max.1ppm(mass/mass)
Zinc (Zn)	max.0.1ppm(mass/mass)
Lead (Pb)	max.0.1ppm(mass/mass)
Iron (Fe)	max.0.2ppm(mass/mass)
Nickel (Ni)	max.0.1ppm(mass/mass)
Substances reducing permanganate (as O)	max.8ppm(mass/mass)

JIS K 8085:2021(\*:Additional test performed by Wako but not required by JIS)The assay of this product is approximately 25%