

SPECIFICATION

Ammonium Nitrate

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|---|
| * Appearance Assay (NH ₄ NO ₃) Solubility in water Residue after ignition (as sulfate) pH of a 50g/l solution at 25°C Chloride (Cl) Nitrite (NO ₂) Phosphate (PO ₄) Sulfate (SO ₄) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe) | White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test max.0.01%(mass/mass) 4.5~6.0 max.3ppm(mass/mass) max.2ppm(mass/mass) max.2ppm(mass/mass) max.0.002%(mass/mass) max.2ppm(mass/mass) max.0.001%(mass/mass) max.3ppm(mass/mass) max.2ppm(mass/mass) |

JIS K 8545:2007(*:Additional test performed by Wako but not required by JIS)