

Code No.01-0333

## SPECIFICATION

## Diammonium Hydrogenphosphate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	7.7~8.3
Chloride (Cl)	max.0.005%
Sulfate (SO <sub>4</sub> )	max.0.01%
Heavy metals (as Pb)	max.0.002%
Arsenic (As)	max.3ppm
Iron (Fe)	max.0.002%
Assay	min.98.5%