

Code No.01-0351

SPECIFICATION

Diammonium (+)-Tartrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Residue after ignition (as sulfate)	max.0.05%
pH of a 50g/L solution at 25°C	6.0~7.0
Chloride (Cl)	max.0.001%
Phosphate (PO ₄)	max.0.001%
Sulfur compounds (as SO ₄)	max.0.005%
Sodium (Na)	max.0.01%
Potassium (K)	max.0.01%
Copper (Cu)	max.1ppm
Zinc (Zn)	max.2ppm
Lead (Pb)	max.1ppm
Iron (Fe)	max.0.001%
Assay	min.99.0%