

Code No.01-0399

## SPECIFICATION

## Aniline

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
Assay(C <sub>6</sub> H <sub>7</sub> N)	min.99.0%(mass/mass)
Density (20°C)	1.019~1.023g/mL
Refractive index 20°C	1.584~1.587
Water	max.0.1%(mass/mass)
Residue after evaporation	max.0.01%(mass/mass)
Nitrobenzene (C <sub>6</sub> H <sub>5</sub> NO <sub>2</sub> ) (GC)	max.0.003%(mass/mass)

JIS K 8042:2022