

Code No.01-0413

## SPECIFICATION

## p-Methoxybenzaldehyde

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellowish brown, clear liquid
Density (20°C)	1.122~1.128g/mL
Refractive index 20°C	1.569~1.574
Water	max.0.5%
Acidity (as CH <sub>3</sub> OC <sub>6</sub> H <sub>4</sub> COOH)	max.1.0%
Assay (capillary GC)	min.97.0%