

Code No.01-0485

SPECIFICATION

Atropine Sulfate Monohydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystalline powder ~ powder
Solubility in water	to pass test
Loss on drying at 110°C	max.5.0%
Residue after ignition	max.0.1%
Acidity (as H ₂ SO ₄)	max.0.03%
Foreign alkaloids	to pass test
Hyoscyamine (after drying)	to pass test
Substances darkened by sulfuric acid	to pass test
Assay (as anhydrous) (after drying)	min.98.0%