

Code No.01-0957

SPECIFICATION

Acetaldehyde

Extra Pure

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
** Assay (CH ₃ CHO)	88.0~92.0%(mass/mass)
Solubility in water	to pass test
Residue after evaporation	max.0.01%(mass/mass)
Acidity (as CH ₃ COOH)	max.0.6%(mass/mass)

JIS K 8030:2010(*:Additional test performed by Wako but not required by JIS)(**:The specification is originally set by Wako not JIS)The assay of this product is approximately 90%