

Code No.01-0957

## SPECIFICATION

## Acetaldehyde

Extra Pure

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
** Assay (CH <sub>3</sub> CHO)	88.0~92.0%(mass/mass)
Solubility in water	to pass test
Residue after evaporation	max.0.01%(mass/mass)
Acidity (as CH <sub>3</sub> COOH)	max.0.6%(mass/mass)

JIS K 8030:2010(\*:Additional test performed by Wako but not required by JIS )(\*\*:The specification is originally set by Wako not JIS)The assay of this product is approximately 90%