

Code No.01-1373

SPECIFICATION

Acetic Anhydride-2,6-Dimethylpyridine-THF Solution

for Nucleic Acid Synthesis

| REQUIREMENT | SPECIFICATION |
|--|--|
| Appearance | Colorless ~ nearly colorless, clear liquid |
| Acetic anhydride (capillary GC) | 8.5~11.5% |
| 2,6-Dimethylpyridine (volumetric analysis) | 9.0~9.6w/v% |
| 2,6-Dimethylpyridine (capillary GC) | 8.0~11.5% |
| Tetrahydrofuran (capillary GC) | 78.0~90.0% |