

Code No.01-1510

## SPECIFICATION

Albumin, from Bovine Serum, Globulin Free  
for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale brown or slightly greenish yellow, crystals ~ crystalline powder to pass test
Solubility in water	
Loss on drying over silica gel in vacuum	max.3.0%
pH of a 10g/L solution at 25°C	6~8
Heavy metals (as Pb)	max.0.002%
Total nitrogen (N)	14~16%
Globulin	no detection
Assay (electrophoresis)	min.98%