

Code No.01-1920

## SPECIFICATION

Acetone for Pesticide Residue And Polychlorinated Biphenyl Tests  
Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance * Suitability for pesticide residue analysis (GC-FPD) Assay (CH <sub>3</sub> COCH <sub>3</sub> ) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis	Colorless clear liquid to pass test min.99.5%(area/area) 0.789~0.792g/mL max.0.3%(mass/mass) max.2ppm(mass/mass) to pass test

JIS K 8040:2020(\*:Additional test performed by Wako but not required by JIS)