

Code No.01-2068

SPECIFICATION

Alisol A

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

| REQUIREMENT | SPECIFICATION |
|--|---------------------------------|
| Description (JP method) | White ~ slightly yellow, powder |
| Optical rotation [α]D20 (JP method) | +86~+106° |
| Purity (JP method) | - |
| Related substances (JP method) | to pass test |
| * Solubility in methanol (after drying) | to pass test |
| * Loss on drying over silica gel | max.2.0% |
| * Assay (HPLC) | min.95.0% |

*:Additional test performed by Wako