

Code No.01-2127

## SPECIFICATION

## Albumin, from Bovine Serum, pH5.2

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White~slightly greenish yellow or pale brown, crystals ~ crystalline powder to pass test
Solubility in water	
Loss on drying over silica gel in vacuum	max.5.0%
pH of a 10g/L solution at 25°C	5.0~5.6
Heavy metals (as Pb)	max.0.005%
Assay (SDS-PAGE)	min.95.0%