

SPECIFICATION

Albumin, from Bovine Serum, Cohn Fraction V, pH7.0 for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	Nearly white ~ pale brown, crystals ~ crystalline powder to pass test
Solubility in water	
Loss on drying over silica gel in vacuum	max.5.0%
pH of a 10g/l solution at 25°C	6.5~7.5
Heavy metals (as Pb)	max.0.005%
Bacterial endotoxins	less than 3EU/mg
Electrophoresis test (SDS-PAGE)	to pass test